

Premier Products for Superior Life Science Research

# NWLSS<sup>TM</sup> Methionine Sulfoxide Reductase A ELISA

Product NWK-MSRAH1 For Research Use Only

ELISA kit for quantification of human methionine sulfoxide reductase A in biological samples.

# **Table of Contents**

Section	Page
Introduction	3
Intended Use	3
Test Principle	3
Specifications	4
Kit Contents	4
Required Materials Not Provided	4
Required Instrumentation	5
Warnings, Precautions, Limitations	5
Storage Instructions	5
Assay Preparation	5
Reagent Preparation	6
Sample Handling/Preparation	7
Standard Curve Preparation	7
Assay Protocol	8
Data Analysis	9
Performance Details	10
References	10
Statement of Limited Warranty	10
Notes	11

#### Introduction:

Methionine Sulfoxide Reductase A (MsrA) plays an essential role in protecting cells against oxidative damage. MsrA reduces methionine sulfoxide (MetO) residues in proteins and free MetO to Methionine (Met). The catalytic activity of MsrA is dependent of bound metal and cofactors but also requires reducing equivalents from either DTT or a thioredoxin-regenerating system. Recent studies indicate that there is a connection between MsrA and Alzheimer's disease in mammals.

#### Intended Use:

The NWLSS™ MsrA ELISA kit is intended to be used for the in vitro quantitative determination of Methionine sulfoxide reductase A (MsrA) in human serum, plasma, cell lysates and cell culture supernatants.

#### Test Principle:

The NWLSS™ Methionine Sulfoxide Reductase A ELISA is a sandwich format Enzyme-Linked Immunosorbent Assay (ELISA). The microtiter plate provided in this kit has been pre-coated with a monoclonal antibody specific to human MsrA. Samples are pipetted into these wells. Nonbound MsrA and other components of the sample are removed by washing, then monoclonal antibody specific to MsrA added. In order to quantitatively determine the amount of MsrA present in the sample, Streptavidin Horseradish Peroxidase (HRP) is added to each microplate well. After another wash step, a TMB-substrate solution added to each well. Finally, a sulfuric acid stop solution is added and the resulting yellow colored product is measured at 450nm. The amount od MsrA in the sample can be determined by direct comparison with the standard curve generated in the assay.

# Specifications:

Format:	1 X 96 well ELISA presented as 12 X 8 well (12 X 8 well) strips in frame.			
Number of tests:	Triplicate = Duplicate =	24 40		
Specificity:	Human MsrA			
Sensitivity:	0.13 ng/mL			
Range:	0.13 ng/mL—5 ng/mL			
<u>Kit Contents:</u> 1 Foil Pouch	96 well microplate precoated with anti-hu MsrA.			
1 vial	Human MsrA Stan	dard (lyophilized)	(1 Vial)	
1 Bottle	Sample/Standard Dilution Buffer		(25mL)	
1 vial	100X Secondary Antibody (Lyophilized) ( (biotin labeled anti-hu MsrA)		(1 Vial)	
1 Bottle	Reagent Dilution Buffer		(25mL)	
1 vial	100X Streptavidin-HRP Conjugate		(150 µL)	
1 Bottle	TMB Substrate Solution (1		(15 mL)	
1 Bottle	Stop Solution (1 N Sulfuric Acid)		(15 mL)	
2 Bottle	20X Concentrated Wash Buffer		(25 mL) ea	
2	Adhesive Plate Co	vers	(2)	

# **Required Materials Not Provided:**

Adjustable micropipettes with disposable tips (5-1000  $\mu L).$  Multi-channel pipettes are useful and help to reduce intra-sample variability.

Polypropylene tubes.

Serological pipettes.

Deionized water.

Automatic plate washer or other aspiration devices are optional.

## **Required Instrumentation:**

Plate reader with **450 nm** capability (650 nm is required for optional monitoring of color development prior to stopping the reaction).

#### Warnings, Precautions & Limitations:

Reagents are intended for research use only and are not for use in diagnostic or therapeutic procedures.

Individual components may be harmful if swallowed, inhaled or absorbed through the skin. Contact should be minimized through the use of gloves and standard good laboratory practices. If contact with skin or eyes occurs, rinse the site immediately with water and consult a physician.

Substrate solutions must be at room temperature prior to use. Avoid contact of substrate solutions with oxidizing agents and metal.

Improper or insufficient washing at any stage of the procedure will result in either false positive or false negative results. Completely empty wells before dispensing fresh Wash Buffer, fill with Wash Buffer as indicated for each wash cycle and do not allow wells to sit uncovered or dry for extended periods.

#### Storage Instructions:

All kit components of this kit are stable at 2 to 8 °C. Any unused reconstituted standard should be discarded or frozen at -70 °C. Standard can be frozen and thawed one time only without loss of immunoreactivity.

#### Assay Preparation:

**1**. Determine the number of wells required to assay standards, samples and controls for the appropriate number of replicates. It is recommended that testing be performed in duplicate or triplicate if possible.

2. Create an assay template showing positioning of standards, controls and samples.

3. Bring all samples and reagents to room temperature before use.

4. To avoid condensation, do not open foil pouches containing the microtiter strips until after they have reached room temperature. Next remove the required number of strips and place in the frame supplied.

Return unused wells to the storage bag with desiccant, seal and store at 2-8  $^{\circ}\mathrm{C}.$ 

# Reagent Preparation:

Secondary Antibody 1. Reconstitute 100X Secondary Antibody by adding 150 µL Reagent Dilution Buffer to the vial.

2. Equilibrate 100X Secondary Antibody to room temperature, mix gently.

3. Mix  $20\mu$ L of 100X Secondary Antibody with 2mL Reagent Dilution Buffer for each 16 well strip to be assayed. Label as **"Working Secondary Antibody Solution".** 

4. Return the unused 100X Secondary Antibody to the refrigerator.

# AVIDIN-HRP Conjugate

1. Equilibrate to room temperature, mix gently.

2. Mix  $20\mu$ L of 100X AVIDIN-HRP Conjugate with 2mL Reagent Dilution Buffer for each16-well strip to be assayed. Label as **"Working Conjugate Solution".** 

3. Return the unused 100X AVIDIN-HRP Conjugate to the refrigerator.

#### Wash Buffer

1. Equilibrate to room temperature, mix to re-dissolve any precipitated salt.

2. Mix 0.5 volume 20X Wash Buffer with 9.5 volumes of deionized water. Label as **"Working Wash Solution".** 

3. Store both the remaining concentrated Wash Buffer and the Working Wash Solution at 4 °C in the refrigerator.

TMB Substrate The TMB Substrate is provided ready to use.

Stop Solution The Stop Solution is provided ready to use

#### Sample Handling/Preparation

The rate of degradation of human MsrA in various matrices has not been fully investigated. It is beyond the scope of this publication to comment on specific sample processing protocols except to recommend that serum or plasma samples be separated by centrifugation and separated from coagulated or packed cells as soon as possible. For best results, serum and plasma samples should be diluted 10X prior to assay.

Cell lysates can be made by isolating cell samples followed by mechanical homogenization or sonication. Homogenates should be clarified by centrifugation prior to assay.

#### Standard Curve Preparation:

Reconstitute the human MsrA standard to 50 ng/mL by adding 1mL of *Standard/Sample Dilution Buffer* into the standard protein glass vial containing lyophilized human MsrA protein. Swirl or mix gently, and allow to sit for 5 minutes to ensure complete reconstitution.

1. Label tubes 1-8 tubes as: 5000, 2500, 1250, 625, 312.5, 156.3, 78.13 and zero (0) ng/mL.

2. Add 900  $\mu L$  Standard Dilution Buffer to tube 1 and 500  $\mu L$  Standard Dilution Buffer to each tube 2-8.

3. Add 100 µL Reconstituted 50 ng/mL Standard to tube 1 and mix well.

4. Make a serial dilution by transferring 500  $\mu$ L of 5000 pg/mL Standard into tube 2 mixing thoroughly then 500  $\mu$ L of resulting 2500 pg/mL to tubes 3 and so on to create all Standards down to 78.13 pg/mL. Sufficient for 1 X 16 well strip (1 set of standards in duplicate).

# Assay Protocol:

1. Add  $100\mu L$  of *Diluted Standards* to the appropriate microtiter wells and 100ul of *Sample Dilution Buffer* to zero wells.

2. Add  $100\mu L$  of **Sample** to each well according to plan.

3. Cover the plate with the plate cover and incubate for 2 hours at 37  $\,^{\circ}\text{C}.$ 

4. Aspirate or decant the solution from the wells then Wash wells 3 times as follows: Dispense 300 uL *Working Wash Solution* to each well and allow to soak for 1-3 minutes before decanting or aspirating the remaining solution from the wells.

5. Add 100µL of *Working Secondary Antibody* to each well.

6. Cover the plate with the plate cover and incubate for 1 hour at 37 °C.

7. Aspirate or decant the solution from the wells then wash the wells 3 times as previously described in step 4.

8. Add 100µL Working Conjugate Solution to each well.

9. Cover the plate with the plate cover and incubate for 30 minutes at room temperature 37  $\,^{\circ}\text{C}.$ 

10. Thoroughly aspirate or decant the solution from the wells. Wash the wells 3 times previously described in step 3.

11. Add 100 $\mu L$  of TMB Substrate to each well. The liquid in the wells should begin to turn blue.

12. Incubate the plate at room temperature for approximately 7-12 minutes.

Note: The incubation time for the TMB substrate is dependent on ambient conditions as well as the specific microtiter plate reader in use. The user should adjust this time as necessary by monitoring the development of blue color at 650 nm and stopping when the high standard has reached maximal absorbance level.

#### Assay Protocol: (continued):

13. After appropriate incubation time, add  $100\mu$ L of **Stop Solution** to each well. The solution in the wells should change from blue to yellow.

14. Read and record the absorbance of each well at 450nm within 20 minutes of adding the Stop Solution.

#### Data Analysis:

1. Plot the mean absorbance at 450 nm for each standard versus the MsrA concentration. Select the best possible fit for the curve obtained (4–parameter is recommended). This can typically be done using the software provided with most plate readers. An example curve is shown below.

2. Sample MsrA is determined by comparing their absorbance at 450 with those of the standard curve.

3. Sample data as read from the standard curve must be multiplied by the dilution factor used.

**Note:** Samples with an ABS<sub>450</sub> exceeding that of the highest standard should be additionally diluted with Sample Dilution Buffer and re-assayed in order to avoid erroneous results.



# Performance Details:

The following substances were tested and found to have no crossreactivity: human MsrB.

## Sensitivity

The minimal detectable dose of human MsrA was calculated to be 0.127 ng/ml, by subtracting three standard deviations from the mean of 12 zero standard replicates (ELISA buffer, S0) and intersecting this value with the standard curve obtained in the same calculation.

Precision Intra-assay = 3.57 % Inter-assay = 3.41 %

Accuracy: Recovery on addition is 100.2~106.6%

(Average 103.3%)

# References:

Davis, D. A. et al. (2000) Biochem. J. 346, 305-311.
Sun, H. et al. (1999) Biochemistry. 38, 105-112.
Gabbita, S. P. et al. (1999) J. Neurochem. 73, 1660-1666.
Abrams, W. R. et al. (1981) Proc. Natl. Acad. Sci. U. S. A. 78, 7483-7486.

# Statement of Limited Warranty:

Northwest Life Science Specialties, LLC (NWLSS) makes no guarantee of any kind, expressed or implied, that extends beyond the description of the material in this kit, except that they will meet our specifications at the time of delivery. Customer's remedy and NWLSS' sole liability is limited to, at NWLSS' option, refund of the purchase price, or the replacement of material not meeting our specification. By acceptance of our product, customer assumes all liability and will indemnify and hold NWLSS harmless for the consequence of this product's use or misuse by the customer, its employees, or others. Refund or replacement is conditioned of customer notifying NWLSS within twenty-one (21) days of the receipt of product. Failure to give notice within 21 days shall constitute a waiver by the customer of all claims hereunder with respect to said product.

# To Order Call: 1-888-449-3091

Notes:



# 5131 NE 94th Avenue, Suite 201 Vancouver, WA 98662 Phone 360-449-3091 or Toll Free: 888-449-3091 Fax 360-449-3092 E-mail: sales@nwlifescience.com