


Northwest
Life Science Specialties, LLC

Anti Myeloperoxidase monoclonal Polyclonal

Catalog Number:	NWA-MPOPA052
Description:	Rabbit IgG anti Hu-Myeloperoxidase (100 µL per Vial)
Specificity:	Specific for human myeloperoxidase
Format:	100 uL IgG in HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Antigen:	Synthetic Peptide
Suggested Use:	Western Blot (1:2000) <i>User should determine optimum titer for each application.</i>
Cross Reactivity:	Human
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment



Immunoblot Analysis of cell lysate
Lane 1 : Purified human MPO protein
Lane 2 : HL-60 cell lysate

- References:**
- 1) Winterbourn CC, et al, Curr Opin Hematol. 200; vol.7(1): pp.53-8.
 - 2) Hansson M, et al, Arch Biochem Biophys. 2006; vol.445(2): pp.214-24.
 - 3) Nicholls SJ, Hazen SL. Arterioscler Thromb Vasc Biol. 2005; vol.25(6): pp.1102-11.
 - 4) Klebanoff SJ. J Leukoc Biol. 2005; vol.77(5): pp.598-625.

For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com