


Northwest
Life Science Specialties, LLC

Anti Peroxiredoxin I Polyclonal

Catalog Number:	NWA-PRX1P086
Description:	Rabbit IgG anti Hu-Peroxiredoxin I (100 µL per Vial)
Specificity:	Specific for human Peroxiredoxin I
Format:	100 uL IgG in HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Antigen:	Recombinant Human Peroxiredoxin I from E. coli.
Suggested Use:	Western Blot (1:2000) <i>User should determine optimum titer for each application.</i>
Cross Reactivity:	Human, Mouse & Rat
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment



IMMUNOBLOT ANALYSIS of cell lysates
Lane 1 : HeLa cell lysate
Lane 2 : A431 cell lysate

- References:**
- (1) Wood, Z. A. et al (2003) Trends Biochem Sci. 28(1):32-40.
 - (2) Angelika Tolle, et al (2005) Free Radic Biol Med. 38:1401-0108
 - (3) Rhee Sue Goo, et al (2001) IUBMB life 52:35-41
 - (4) Min Hee Choi, et al (2005) Nature letters 435(19):347-353

For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com