


Northwest

Life Science Specialties, LLC

Anti EC Superoxide Dismutase (Extracellular: SOD3) Monoclonal Clone 1H12

Catalog Number:	NWA-SOD31H12
Description:	Murine IgG anti hu-Superoxide Dismutase 3 (100 µL per Vial).
Specificity:	Specific for Human SOD3
Format:	100 µL IgG _{2bκ} in HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Antigen:	Recombinant human protein purified from E.coli
Reconstitution:	NA
Suggested Use:	ELISA Western Blot (1:2000) Immunoprecipitation (1ul) Immunohistochemistry User should determine optimum titer for each Application.
Cross Reactivity:	Human
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment
References:	1) Kuninaka, S. et al. (2000) Br. J. Cancer. 83, 928-934. 2) Strange, R. W. et al. (2003) J. Mol. Biol. 328,877-891. 3) Weisiger, R. A., and Fridovich, I. (1973) J. Biol. chem. 248, 3582-3592. 4) Enghild, J. J. et al. (1999) Biochem J. 317, 51-57.



For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com